

## Datasheet: AHP2281

<b>Description:</b>	RABBIT ANTI Bcl-2 INTERACTING KILLER
<b>Specificity:</b>	Bcl-2 INTERACTING KILLER
<b>Other names:</b>	BLK
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			0.5 - 4.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified IgG - liquid
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	50% Glycerol
<b>Stabilisers</b>	1% Bovine Serum Albumin 0.02% Thiomersal
<b>Approx. Protein Concentrations</b>	IgG concentration 0.2 mg/ml
<b>Immunogen</b>	Synthetic peptide surrounding amino acid 30 of mouse Bcl-2 interacting killer
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">O70337</a> <a href="#">Related reagents</a>

<b>Specificity</b>	<b>Rabbit anti Bcl-2 interacting killer antibody</b> recognizes Bcl-2-interacting killer, also known as BLK. Bcl-2 interacting killer plays a role in mediating apoptotic responses via interaction with Bcl-2 and Bcl-xL.
<b>Western Blotting</b>	Rabbit anti Bcl-2 Interacting kiler antibody recognizes a band of approximately 27 kDa in mouse kidney tissue lysate
<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10585 available at: 10585: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10585.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10585.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M345558:190125'

Printed on 11 Jul 2021