

Datasheet: AHP2301

Description:	RABBIT ANTI CYR61
Specificity:	CYR61
Other names:	INSULIN-LIKE GROWTH FACTOR-BINDING PROTEIN 10
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			5 ug/ml
Western Blotting	▪			0.5 - 4.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Species Cross Reactivity

Reacts with: Human

Based on sequence similarity, is expected to react with:Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

30% Glycerol

0.5% Bovine Serum Albumin

0.01% Thiomersal

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Synthetic peptide surrounding amino acid 166 of mouse CYR61
External Database Links	<p>UniProt: P18406 Related reagents</p> <p>Entrez Gene: 16007 Cyr61 Related reagents</p>
Synonyms	Ccn1, Igfbp10
Specificity	Rabbit anti CYR61 antibody recognizes cysteine-rich angiogenic inducer 61 (CYR61) also known as insulin-like growth factor-binding protein 10 (IGFBP-10/ IBP-10). CYR61 interacts with integrin receptors thereby regulating a variety of processes including cell proliferation and cell adhesion. CYR61 has also been shown to play a role in angiogenesis, inflammation and chemotaxis.
Western Blotting	Rabbit anti CYR61 antibody recognizes a band of approximately 42 kDa in Jurkat cell lysate and mouse small intestine tissue lysates
Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10096 available at: https://www.bio-rad-antibodies.com/SDS/AHP2301 10096
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M429926:240501'

