

## Datasheet: AHP2316

<b>Description:</b>	RABBIT ANTI HEME OXYGENASE 2
<b>Specificity:</b>	HEME OXYGENASE 2
<b>Other names:</b>	HO-2 .
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			20 ug/ml
Immunoprecipitation	▪			20 ug/ml
Western Blotting	▪			0.5 - 4.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Rat

Based on sequence similarity, is expected to react with: Mouse, Monkey, Hamster  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

30% Glycerol  
 0.5% Bovine Serum Albumin  
 0.01% Thiomersal

<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Synthetic peptide surrounding amino acid 34 of human heme oxygenase 2
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P30519</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">3163</a>    HMOX2    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	HO2
<b>Specificity</b>	<b>Rabbit anti heme oxygenase 2 antibody</b> recognizes heme oxygenase 2 (HO-2). Heme oxygenases cleaves heme into the bilirubin precursor biliverdin. HO-2 represents the constitutive form of the enzyme.
<b>Western Blotting</b>	Rabbit anti heme oxygenase 2 antibody recognizes a band of approximately 36kDa in rat kidney tissue lysate
<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10096 available at: 10096: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10096.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10096.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...)    [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...)    [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...)    [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M337254:181217'

