

Datasheet: AHP2376

BATCH NUMBER 153221

Description:	RABBIT ANTI BOVINE INTERLEUKIN-17A
Specificity:	IL-17A
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1.0 ug/ml - 10 ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 ug/ml - 1.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Dolphin, Mouse, Sheep, Pig
Reacts weakly with: Dog, Horse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antisera to bovine interleukin-17A were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine interleukin-17A
External Database Links	<p>UniProt: Q687Y7 Related reagents</p> <p>Entrez Gene: 282863 IL17A Related reagents</p>
Specificity	<p>Rabbit anti Bovine interleukin-17A antibody recognizes bovine interleukin-17A (IL-17A), a 35 kDa disulfide-linked homodimer consisting of two 153 amino acid chains, which belongs to the IL-17 family (classified as IL-17A to IL-17F), and signals through the IL-17 receptor (IL-17R/CDw217).</p> <p>IL-17, originally thought to have been produced primarily by activated CD4+ T cells, appears to have a range of cellular origins and induces a variety of target cells to secrete inflammatory cytokines, including IL-1β, IL-6, TNFα, IFNγ and granulocyte colony-stimulating factor, as well as synergizing with other inflammatory cytokines and agonists.</p>
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with AHP2376B as the detection reagent. Recombinant bovine interleukin-17A (PBP025) may be used as a standard.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP237610040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[RECOMBINANT BOVINE INTERLEUKIN-17A \(PBP025\)](#)

[RABBIT ANTI BOVINE INTERLEUKIN-17A:Biotin \(AHP2376B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M360835:200122'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)