

Datasheet: AHP2378B

Description:	GOAT ANTI BOVINE INTERLEUKIN-1RA:Biotin
Specificity:	IL-1RA
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.01ug/ml - 0.5ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to bovine interleukin-1RA were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine interleukin-1RA
External Database Links	UniProt: O77482 Related reagents

Entrez Gene:[281860](#) [IL1RN](#) [Related reagents](#)

Specificity	Goat anti Bovine interleukin-1RA antibody recognizes bovine interleukin-1 receptor antagonist, a 19.93 kDa protein belonging to the IL-1 cytokine family. Bovine IL-1RA is secreted by several cell types, including epithelial cells and adipocytes. IL-1RA functions by inhibiting the pro-inflammatory effects of both IL-1 α and IL-1 β .
ELISA	This product may be used as a detection antibody in a sandwich ELISA together with AHP2378 as the capture reagent. Recombinant Bovine interleukin-1RA (PBP033) may be used as a standard.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M334175:181127'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)