

Datasheet: AHP2382

Description:	RABBIT ANTI BOVINE VEGF-A
Specificity:	VEGF-A
Other names:	VPF
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1.0 ug/ml - 10 ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			1.0 ug/ml - 10 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Species Cross Reactivity	<p>Reacts with: Cat, Dog, Horse, Mouse, Sheep, Pig, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to bovine VEGF-A were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine VEGF-A
External Database Links	<p>UniProt: P15691 Related reagents</p> <p>Entrez Gene: 281572 VEGFA Related reagents</p>
Synonyms	VEGF
Specificity	<p>Rabbit anti Bovine VEGF-A antibody recognizes bovine vascular endothelial growth factor A (VEGF-A), also known as vascular permeability factor (VPF). VEGF-A is a 22.3 kDa disulfide-linked homodimer belonging to the PDGF / VEGF family of growth factors.</p> <p>Bovine VEGF-A is a growth factor that induces cell proliferation in endothelial cells. It is important in endothelial cell growth, vasculogenesis and angiogenesis. VEGF-A is also active in the induction of blood vessel permeability.</p> <p>VEGF-A is well conserved amongst mammals and Rabbit anti Bovine VEGF-A antibody recognizes VEGF-A in a range of species.</p>
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with AHP2382B as the detection reagent. Recombinant bovine VEGF-A (PBP027) may be used as a standard.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP2382 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Recommended Useful Reagents

[RECOMBINANT BOVINE VEGF-A \(PBP027\)](#)
[RABBIT ANTI BOVINE VEGF-A:Biotin \(AHP2382B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M334988:181203'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)