

Datasheet: AHP2392B

Description:	RABBIT ANTI PIG INTERLEUKIN-8:Biotin
Specificity:	IL-8
Other names:	CXCL-8
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.01ug/ml - 0.5ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to pig interleukin-8 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen Recombinant pig interleukin-8

External Database Links

UniProt:
[P26894](#) [Related reagents](#)

Entrez Gene:
[396880](#) IL8 [Related reagents](#)

Synonyms CXCL8

Specificity **Rabbit anti Pig interleukin-8 antibody** recognizes porcine interleukin-8 (IL-8), also known as CXCL-8. Porcine IL-8 is an 11.63 kDa cytokine belonging to the chemokine C-X-C family.

Porcine IL-8 is a chemotactic cytokine secreted by a range of cell types such as macrophages, endothelial cells and epithelial cells in response to inflammation. IL-8 induces chemotaxis of neutrophils and granulocytes to sites of infection.

ELISA This product may be used as a detection antibody in a sandwich ELISA together with [AHP2392](#) as the capture reagent. Recombinant Pig interleukin-8 ([PBP024](#)) may be used as a standard.

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2392B>
10040

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M361221:200210'

Printed on 18 Jan 2024
