

Datasheet: AHP2413

Description:	RABBIT ANTI STAT5A
Specificity:	STAT5A
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Mouse

Species Cross Reactivity

Reacts with: Human, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Purified IgG prepared by affinity chromatography

Antiserum Preparation

Antiserum to STAT5A was raised by repeated immunisation of rabbits with highly purified antigen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.1% Sodium Azide (NaN₃)
0.2% Gelatin

Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Peptide mapping to the carboxy terminus of mouse STAT5A
External Database Links	<p>UniProt: P42230 Related reagents</p> <p>Entrez Gene: 20850 Stat5a Related reagents</p>
Synonyms	Mgf, Mpf
Specificity	<p>Rabbit anti STAT5A antibody recognizes signal transducer and activator of transcription 5A (STAT5A), a member of the STAT protein family. STAT5A plays an important role in cell growth, survival, differentiation and transformation by transducing signals from cytokine receptors to the nucleus.</p> <p>STAT5A is tyrosine phosphorylated by members of the Janus kinase (JAK) family in response to cytokines and growth factor stimulation. Phosphorylation results in STAT activation and the formation of STAT5 homo- and STAT5/STAT3 heterodimers. STAT dimers translocate from the cytoplasm to the nucleus where they act as transcription factors. Deregulated JAK/STAT signaling has been implicated in a variety of cancers, and STATs can regulate the tumor microenvironment to trigger a shift from cancer inhibition to promotion.</p>
Western Blotting	Rabbit anti STAT5A recognizes a band of approximately 91 kDa when run under reducing conditions
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20341 available at: 20341: https://www.bio-rad-antibodies.com/uploads/MSDS/20341.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M361229:200210'

Printed on 23 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)