

Datasheet: AHP2419

Description:	RABBIT ANTI CASPASE-9 p10
Specificity:	CASPASE-9 p10
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Mouse, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography
Antiserum Preparation	Antiserum to caspase-9 p10 was raised by repeated immunisation of rabbits with highly purified antigen
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	<p>0.1% Sodium Azide (NaN₃)</p> <p>0.2% Gelatin</p>

Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Peptide mapping to an internal region of human caspase-9 p10 subunit
External Database Links	<p>UniProt: P55211 Related reagents</p> <p>Entrez Gene: 842 CASP9 Related reagents</p>
Synonyms	MCH6
Specificity	<p>Rabbit anti caspase-9 p10 antibody recognizes caspase-9, a member of the cysteine-aspartic acid protease family. The active form of caspase-9 is generated by cleavage of the pro-caspase-9 protein into fragments, known as the large and small fragments. The trigger for caspase-9 cleavage is the formation of the apoptosome (comprising of cytochrome c, APAF-1, dATP and pro-caspase-9), upon stimulation of the intrinsic apoptosis pathway.</p> <p>The rabbit anti caspase-9 antibody was raised against an internal peptide that is present in one of the small fragments (p10) generated in the pro-caspase-9 cleavage process at Asp330 (Zou et al. 2002). The antibody will therefore detect pro-caspase-9 (as the motif is present) and the p10 caspase-9 fragment. The antibody can therefore be used to distinguish the p10 caspase-9 fragment from other pro-caspase-9 cleavage products.</p>
Western Blotting	Rabbit anti caspase-9 p10 recognizes a band of approximately 47 kDa when run under reducing conditions without induction of apoptosis
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20341 available at: 20341: https://www.bio-rad-antibodies.com/uploads/MSDS/20341.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M345565:190125'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)