

Datasheet: AHP2421

Description:	RABBIT ANTI BAX
Specificity:	BAX
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/400

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Purified IgG prepared by affinity chromatography

Antiserum Preparation

Antiserum to BAX was raised by repeated immunisation of rabbits with highly purified antigen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.1% Sodium Azide (NaN₃)
0.2% Gelatin

Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Peptide mapping to the amino terminus of human BAX
External Database Links	<p>UniProt: Q07812 Related reagents</p> <p>Entrez Gene: 581 BAX Related reagents</p>
Synonyms	BCL2L4
Specificity	<p>Rabbit anti BAX antibody recognizes the Bcl-2-associated X protein (BAX) also known as Bcl-2 like protein 4. BAX together with BAK (BCL-2 antagonist or killer), another member of the Bcl-2 family, plays a key role in the intrinsic apoptosis pathway by inducing mitochondrial outer membrane permeabilization (MOMP) (Giansanti and Scovassi 2008). This event is aided by BID (BH3-interacting domain death agonist), which activates the two proteins and thereby facilitates the formation of homo- or heterodimers. Often in the literature, MOMP is regarded as the point of no return upon which a cell is committed to undergo apoptosis (Chipuk et al. 2006).</p>
Western Blotting	Rabbit anti BAX recognizes a band of approximately 23 kDa when run under reducing conditions
References	1. Olejnik, J. <i>et al.</i> (2013) Ebola virus does not block apoptotic signaling pathways. J Virol. 87 (10): 5384-96.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #20341 available at: https://www.bio-rad-antibodies.com/SDS/AHP2421</p> <p>20341</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M376901:210212'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)