

Datasheet: AHP2432

**BATCH NUMBER 169331**

|                      |                     |
|----------------------|---------------------|
| <b>Description:</b>  | RABBIT ANTI AP2A1   |
| <b>Specificity:</b>  | AP2A1               |
| <b>Format:</b>       | Purified            |
| <b>Product Type:</b> | Polyclonal Antibody |
| <b>Isotype:</b>      | Polyclonal IgG      |
| <b>Quantity:</b>     | 50 µl               |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                  | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪   |    |                | 1/500 - 1/2000     |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antiserum to AP2A1 was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.02% Sodium Azide  
50% Glycerol

**Immunogen** Recombinant human AP2A1

---

**External Database**

**Links**

**UniProt:**

[O95782](#) [Related reagents](#)

**Entrez Gene:**

[160](#) AP2A1 [Related reagents](#)

---

**Synonyms**

ADTAA, CLAPA1

---

**Specificity**

**Rabbit anti AP2A1 antibody** recognizes AP-2 complex subunit alpha-1 (encoded by AP2A1), also known as adaptor related protein complex 2 alpha 1 subunit, clathrin assembly protein complex 2 alpha-A large chain or alpha1-adaptin.

AP-2 complex subunit alpha-1 is a subunit of the adaptor protein complex (AP-2) and is involved in protein transport.

---

**Western Blotting**

Rabbit anti AP2A1 antibody detects a band of approximately 108 kDa in cell lysates and tissue lysates under reducing conditions

---

**Storage**

Store at -20°C

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10049 available at:  
<https://www.bio-rad-antibodies.com/SDS/AHP2432>  
10049

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M354169:190624'

**Printed on 23 May 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)