

## Datasheet: AHP2449

<b>Description:</b>	RABBIT ANTI CD126/IL-6R
<b>Specificity:</b>	CD126/IL-6R
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
-----------------------	-------

<b>Product Form</b>	Purified IgG - liquid
---------------------	-----------------------

<b>Antiserum Preparation</b>	Antiserum to CD126/IL-6R was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
------------------------------	--

<b>Buffer Solution</b>	Phosphate buffered saline
------------------------	---------------------------

<b>Preservative</b>	0.05% ProClin 300
<b>Stabilisers</b>	50% Glycerol

<b>Immunogen</b>	Recombinant human CD126/IL-6R
------------------	-------------------------------

### External Database Links

#### UniProt:

[P08887](#)

[Related reagents](#)

#### Entrez Gene:

[3570](#)

IL6R

[Related reagents](#)

<b>Specificity</b>	<b>Rabbit anti CD126/IL-6R antibody</b> recognizes CD126/IL-6R, also known as interleukin-6 receptor subunit alpha. CD126/IL-6R is a member of the type I cytokine receptor family. CD126/IL-6R binds to IL-6, the IL-6/IL-6R complex then associates with two molecules of interleukin 6 signal transducer (IL6ST), leading to signaling via the JAK/STAT pathway.
<b>Western Blotting</b>	Rabbit anti CD126/IL-6R antibody detects a band of approximately 52 kDa in cell lysates under reducing conditions
<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	ProClin is a trademark of The Dow Chemical Company (“Dow”) or an affiliated company of Dow.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20376 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP2449">https://www.bio-rad-antibodies.com/SDS/AHP2449</a> 20376
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M432705:240912'

Printed on 12 Sep 2024