

## Datasheet: AHP2459

**BATCH NUMBER 141419**

<b>Description:</b>	RABBIT ANTI DFF-45
<b>Specificity:</b>	DFF-45
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			1/50 - 1/200
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Target Species

Human

#### Species Cross Reactivity

Based on sequence similarity, is expected to react with: Mouse, Rat  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid

#### Antiserum Preparation

Antiserum to DFF-45 was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline

#### Preservative Stabilisers

0.02% Sodium Azide  
 50% Glycerol

<b>Immunogen</b>	Recombinant human DFF-45
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">O00273</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">1676</a>    DFFA    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	DFF1, DFF45
<b>Specificity</b>	<p><b>Rabbit anti DFF-45 antibody</b> recognizes DFF, also known as DNA fragmentation factor subunit alpha, DFF1, DFFA, inhibitor of caspase-activated DNase or ICAD.</p> <p>DNA fragmentation factor (DFF) is a heterodimeric protein of 40kDa (DFF-40) and 45 kDa (DFF-45) subunits. During apoptosis caspase-3 cleaves DFF into the DFF-45 and DFF-40 fragments.</p>
<b>Western Blotting</b>	Rabbit anti DFF-45 antibody detects a band of approximately 43 kDa in cell lysates and tissue lysates under reducing conditions
<b>Storage</b>	<p>Store at -20°C</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10049 available at:  <a href="https://www.bio-rad-antibodies.com/SDS/AHP2459">https://www.bio-rad-antibodies.com/SDS/AHP2459</a>  10049</p>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...)    [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...)    [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...)    [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

'M336799:181214'

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)