

Datasheet: AHP2461

BATCH NUMBER 141421

Description:	RABBIT ANTI DNAJA3
Specificity:	DNAJA3
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Mouse, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to DNAJA3 was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	<p>0.02% Sodium Azide</p> <p>50% Glycerol</p>

Immunogen Recombinant human DNAJA3

External Database

Links

UniProt:

[Q96EY1](#)

[Related reagents](#)

Entrez Gene:

[9093](#)

DNAJA3

[Related reagents](#)

Synonyms

HCA57, TID1

Specificity

Rabbit anti DNAJA3 antibody recognizes DNAJA3, also known as DnaJ homolog subfamily A member 3 or tumorous imaginal disc 1(TID1). DNAJA3 is a member of the DNAJ/Hsp40 protein family and as such contains a conserved DnaJ domain. Members of the DNAJ/Hsp40 family interact with and facilitate the activity of the Hsp70 family of molecular chaperones.

Western Blotting

Rabbit anti DNAJA3 antibody detects a band of approximately 50 kDa in cell lysates and tissue lysates under reducing conditions

Storage

Store at -20°C

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10049 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP2461>
10049

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)

[FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...)

[HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...)

[Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...)

[RPE](#)

Sheep Anti Rabbit IgG (STAR36...)

[DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M336814:181214'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)