

Datasheet: AHP2542

**BATCH NUMBER 162385**

<b>Description:</b>	RABBIT ANTI 4EBP1 (pSer65)
<b>Specificity:</b>	4EBP1 (pSer65)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Human

### Species Cross Reactivity

Based on sequence similarity, is expected to react with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antiserum to 4EBP1 (pSer65) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.02% Sodium Azide  
50% Glycerol

**Immunogen** Phospho specific-peptide corresponding to residues surrounding serine 65 of human 4EBP1

**External Database Links**

**UniProt:**

[Q13541](#) [Related reagents](#)

**Entrez Gene:**

[1978](#) EIF4EBP1 [Related reagents](#)

**Specificity**

**Rabbit anti 4EBP1 (pSer65) antibody** recognizes eukaryotic translation initiation factor 4E-binding protein 1 (4EBP1) also known as phosphorylated heat-and acid-stable protein regulated by insulin 1 (PHAS-I), when phosphorylated at serine 65. 4EBP1 belongs to the eIF4E-binding protein family.

4EBP1, when bound to eIF4E, inhibits cap-dependent translation. Phosphorylation at serine 65 inhibits binding of 4EBP1 to eIF4F allowing cap-dependent translation.

**Storage**

Store at -20°C

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

**Guarantee**

12 months from date of despatch

**Health And Safety Information**

Material Safety Datasheet documentation #10049 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2542>  
10049

**Regulatory**

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M376917:210212'

