

Datasheet: AHP2546

Description:	RABBIT ANTI AKT1 (pThr450)
Specificity:	AKT1 (pThr450)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/500 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Human

Product Form Purified IgG - liquid

Antiserum Preparation Antiserum to AKT1 (pThr450) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative 0.02% Sodium Azide
Stabilisers 50% Glycerol

Immunogen Phospho specific-peptide corresponding to residues surrounding threonine 450 of human AKT1

External Database Links

UniProt:
[P31749](#) [Related reagents](#)

Entrez Gene:
[207](#) AKT1 [Related reagents](#)

Synonyms PKB, RAC

Specificity **Rabbit anti AKT1 (pThr450) antibody** recognizes AKT1, also known as RAC-alpha serine/threonine-protein kinase, protein kinase B alpha or Proto-oncogene c-Akt, when phosphorylated at threonine 450. AKT1 is a 480 amino acid ~48 kDa protein kinase that play a key

role in the PI3K-Akt signaling pathway, known to regulate cell growth, proliferation and survival. Phosphorylation of threonine 308 and serine 473 are required for full activation of AKT1. Newly synthesized Akt is suggested to require phosphorylation at threonine 450 for proper protein folding ([Liao and Hung 2010](#)).

Further Reading	1. Liao Y, Hung M-C. (2010) Physiological regulation of Akt activity and stability American Journal of Translational Research. 2 (1): 19-42.
------------------------	--

Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10049 available at: 10049: https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M360493:191126'

Printed on 01 May 2020