

Datasheet: AHP2564

BATCH NUMBER 141526

Description:	RABBIT ANTI CATENIN BETA (pTyr333)
Specificity:	CATENIN BETA (pTyr333)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Human

Product Form Purified IgG - liquid

Antiserum Preparation Antiserum to catenin beta (pTyr333) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution Phosphate buffered saline

Preservative 0.02% Sodium Azide
Stabilisers 50% Glycerol

Immunogen Phospho specific-peptide corresponding to residues surrounding tyrosine 333 of human catenin beta

External Database Links

UniProt:
[P35222](https://www.uniprot.org/entry/P35222) [Related reagents](#)

Entrez Gene:

[1499](#) CTNNB1 [Related reagents](#)

Synonyms	CTNNB
-----------------	-------

Specificity	<p>Rabbit anti catenin beta (pTyr333) antibody recognizes catenin beta, when phosphorylated at tyrosine 333. Catenin beta is a component of adherens junctions protein complexes that are involved in cell adhesion and actin cytoskeleton organization. Catenin beta also plays an important role in the canonical Wnt signaling pathway.</p> <p>Protein tyrosine kinase 6 (PTK6) phosphorylates catenin beta at tyrosine 64, tyrosine 142, tyrosine 331 and/or tyrosine 333.</p>
--------------------	---

Further Reading	1. Palka-Hamblin, H.L. <i>et al.</i> (2010) Identification of beta-catenin as a target of the intracellular tyrosine kinase PTK6 J Cell Sci. 123 (Pt 2): 236-45.
------------------------	--

Storage	<p>Store at -20°C</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	---

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10049 available at: https://www.bio-rad-antibodies.com/SDS/AHP2564 10049
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@680](#), [DyLight@800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M361237:200210'

