

Datasheet: AHP2566

Description:	RABBIT ANTI CD117 (pTyr721)
Specificity:	CD117 (pTyr721)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Based on sequence similarity, is expected to react with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antiserum to CD117 (pTyr721) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.02% Sodium Azide
50% Glycerol

Immunogen

Phospho specific-peptide corresponding to residues surrounding tyrosine 721 of human CD117

External Database**Links****UniProt:**

[P10721](#) [Related reagents](#)

Entrez Gene:

[3815](#) KIT [Related reagents](#)

Specificity

Rabbit anti CD117 (pTyr721) antibody recognizes CD117, also known as proto-oncogene c-kit, mast/stem cell growth factor receptor Kit, Piebald trait protein and v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog, when phosphorylated at tyrosine 721. CD117 is a member of the tyrosine protein kinase family and is the receptor for stem cell factor (SCF).

Upon binding of SCF, CD117 activates a number of cell signaling pathways including MAP kinase, PI3K and JAK/STAT pathways ([Laing *et al.* 2013](#)).

SCF binding induces autophosphorylation of CD117 at tyrosine 568, tyrosine 570, tyrosine 721 and tyrosine 730. Phosphorylation of tyrosine 721 of CD117 is important for the binding and activation of PI-3 kinase p85 subunit alpha ([Lennartsson and Rönstrand, 2012](#)).

Further Reading

1. Liang, J. *et al.* (2013) The C-kit receptor-mediated signal transduction and tumor-related diseases [Int J Biol Sci. 9 \(5\): 435-43](#).
 2. Lennartsson, J. & Rönstrand, L. (2012) Stem cell factor receptor/c-Kit: from basic science to clinical implications [Physiol Rev. 92 \(4\): 1619-49](#).
-

Storage

Store at -20°C

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10049 available at: 10049: <https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

TidyBlot WESTERN BLOT DETECTION REAGENT:HRP (STAR209P)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M344130:190111'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)