

## Datasheet: AHP2583

<b>Description:</b>	RABBIT ANTI C-MYC (pThr358)
<b>Specificity:</b>	C-MYC (pThr358)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antiserum to c-Myc (pThr358) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.02% Sodium Azide
<b>Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Phospho specific-peptide corresponding to residues surrounding threonine 358 of human c-Myc
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P01106</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b>

---

<b>Synonyms</b>	BHLHE39
-----------------	---------

---

<b>Specificity</b>	<p><b>Rabbit anti c-Myc (pThr358) antibody</b> recognizes proto-oncogene c-Myc, also known as Class E basic helix-loop-helix protein 39, transcription factor p64 and avian myelocytomatosis viral oncogene homolog, when phosphorylated at threonine 358. c-Myc is a transcription factor and a member of the Myc family of basic helix-loop-helix-leucine zipper (bHLH-Zip) transcription factors. c-Myc is involved in regulation of the cell cycle, cell differentiation and apoptosis.</p> <p>c-Myc is phosphorylated at threonine 358 by PAK2, this negatively regulates c-Myc activity by inhibiting DNA binding (<a href="#">Huang et al. 2004</a>).</p>
--------------------	--

---

<b>Further Reading</b>	1. Huang, Z. <i>et al.</i> (2004) Negative control of the Myc protein by the stress-responsive kinase Pak2. <a href="#">Mol Cell Biol. 24 (4): 1582-94.</a>
------------------------	---

---

<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10049 available at: 10049: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M336640:181214'

