

Datasheet: AHP2594

BATCH NUMBER 141554

Description:	RABBIT ANTI EGF R (pTyr1148)
Specificity:	EGF R (pTyr1148)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to EGF R (pTyr1148) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	0.02% Sodium Azide
Stabilisers	50% Glycerol
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 1148 of human EGF R

External Database Links

UniProt:
[P00533](https://www.uniprot.org/entry/P00533) [Related reagents](#)

Entrez Gene:[1956](#) EGFR [Related reagents](#)**Synonyms**

ERBB1

Specificity

Rabbit anti EGF R (pTyr1148) antibody recognizes the epidermal growth factor receptor (EGF R), also known as proto-oncogene c-ErbB-1 or receptor tyrosine-protein kinase erbB-1, when phosphorylated at tyrosine 1148. EGF R is a transmembrane glycoprotein and a member of the ErbB receptor tyrosine kinase family. EGF R is a cell surface receptor for members of the epidermal growth factor family.

EGF R plays a role in numerous cellular processes such as cell proliferation and apoptosis through the activation of downstream signaling pathways ([Scaltriti and Baselga 2006](#)).

EGF R is autophosphorylated at tyrosine 1148 in response to ligand binding and receptor dimerization.

Further Reading

1. Scaltriti, M. & Baselga, J. (2006) The epidermal growth factor receptor pathway: a model for targeted therapy. [Clin Cancer Res. 12 \(18\): 5268-72.](#)

Storage

Store at -20°C

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10049 available at:
10049: <https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)Sheep Anti Rabbit IgG (STAR35...) [RPE](#)Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M360523:191126'

Printed on 24 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)