

## Datasheet: AHP2599

<b>Description:</b>	RABBIT ANTI ELK-1 (pSer383)
<b>Specificity:</b>	ELK-1 (pSer383)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			1/50 - 1/100
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Based on sequence similarity, is expected to react with: Mouse, Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species.
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antiserum to ELK-1 (pSer383) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.02% Sodium Azide
<b>Stabilisers</b>	50% Glycerol
<b>Immunogen</b>	Phospho specific-peptide corresponding to residues surrounding serine 383 of human ELK-1
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P19419</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">2002</a>   ELK1   <a href="#">Related reagents</a></p>
<b>Specificity</b>	<b>Rabbit anti ELK-1 (pSer383) antibody</b> recognizes ELK-1, also known as ETS domain-containing

protein Elk-1, when phosphorylated at serine 383.

ELK-1 is a member of the ternary complex factor (TCF) family of transcription factors, and plays a role in mediating gene expression in response to growth factors.

ELK-1 is phosphorylated at serine 383 and 389 in its C-terminal region by mitogen activated protein (MAP) kinase. Phosphorylation of serine 383 and 389 is important for activation of ELK1 ([Li et al. 2003](#)).

---

<b>Further Reading</b>	1. Li, Q.J. <i>et al.</i> (2003) MAP kinase phosphorylation-dependent activation of Elk-1 leads to activation of the co-activator p300. <a href="#">EMBO J. 22 (2): 281-91.</a>
<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10049 available at: 10049: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf</a>
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)  
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)  
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)  
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)  
Sheep Anti Rabbit IgG (STAR36...) [DyLight@488](#), [DyLight@549](#), [DyLight@649](#),  
[DyLight@680](#), [DyLight@800](#)  
Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

'M336911:181214'

Printed on 16 Mar 2019