

Datasheet: AHP2688

Description:	RABBIT ANTI STAT5A (pSer780)
Specificity:	STAT5A (pSer780)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antiserum to STAT5A (pSer780) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative	0.02% Sodium Azide
Stabilisers	50% Glycerol

Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 780 of human STAT5A
------------------	---

External Database Links

UniProt:

[P42229](https://www.uniprot.org/uniprot/P42229)

[Related reagents](#)

Entrez Gene:[6776](#) STAT5A [Related reagents](#)

Synonyms	STAT5
-----------------	-------

Specificity	Rabbit anti STAT5A (pSer780) antibody recognizes STAT5A, also known as when phosphorylated at serine 780. STAT5A is a member of the STAT family of transcription factors. Members of the STAT family reside in the cytosol until activated. Upon activation STATs dimerize and translocate to the nucleus where they exert their transcriptional effects. Serine residue 780 is located in the transcriptional activation domain of STAT5A.
--------------------	--

Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10049 available at: https://www.bio-rad-antibodies.com/SDS/AHP2688 10049
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)Sheep Anti Rabbit IgG (STAR35...) [RPE](#)Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M361231:200210'

Printed on 12 Aug 2023