

Datasheet: AHP2708

Description:	RABBIT ANTI CATHEPSIN H
Specificity:	CATHEPSIN H
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
-----------------------	-------

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antiserum to cathepsin H was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative	0.02% Sodium Azide (NaN ₃)
Stabilisers	0.5% Bovine Serum Albumin 50% Glycerol

Immunogen	Synthesized peptide derived from the internal region of human cathepsin H (amino acids 70-150)
------------------	--

External Database Links

UniProt:

[P09668](https://www.uniprot.org/uniprot/P09668)

[Related reagents](#)

Entrez Gene:

[1512](#) CTSH [Related reagents](#)

Synonyms	CPSB
-----------------	------

Specificity	Rabbit anti Cathepsin H antibody recognizes pro-cathepsin H, also known as CTSH and CPSB. Cathepsin H is a lysosomal cysteine protease that is synthesized as an inactive proenzyme composed of prodomain and catalytic domains. It is currently unclear if the processing/cleavage of the proenzyme is performed by an autocatalytic method (Rojnik et al. 2012) or requires other enzymes for its activation (Hao et al. 2018). Cathepsin H has been reported to be overexpressed in human prostate tumor cells (Waghray et al. 2002).
--------------------	---

Western Blotting	Rabbit anti cathepsin H antibody detects a band of approximately 40 kDa in HEK 293 cell lysate
-------------------------	--

Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10048 available at: 10048: https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M353137:190410'

Printed on 02 Mar 2021