

## Datasheet: AHP2716

**BATCH NUMBER 152291**

<b>Description:</b>	RABBIT ANTI BAX
<b>Specificity:</b>	BAX
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
-----------------------	-------

<b>Product Form</b>	Purified IgG - liquid
---------------------	-----------------------

**Antiserum Preparation** Antiserum to BAX was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

<b>Buffer Solution</b>	Phosphate buffered saline
------------------------	---------------------------

<b>Preservative</b>	0.02% Sodium Azide
<b>Stabilisers</b>	50% Glycerol

<b>Immunogen</b>	Synthetic human BAX peptide
------------------	-----------------------------

#### External Database Links

##### UniProt:

[Q07812](#)

[Related reagents](#)

##### Entrez Gene:

---

<b>Synonyms</b>	BCL2L4
-----------------	--------

---

<b>Specificity</b>	<b>Rabbit anti BAX antibody</b> recognizes the Bcl-2-associated X protein (BAX) also known as Bcl-2 like protein 4. BAX together with BAK (BCL-2 antagonist or killer), another member of the Bcl-2 family, plays a key role in the intrinsic apoptosis pathway by inducing mitochondrial outer membrane permeabilization (MOMP) ( <a href="#">Giansanti and Scovassi 2008</a> ). This event is aided by BID (BH3-interacting domain death agonist), which activates the two proteins and thereby facilitates the formation of homo- or heterodimers. Often in the literature, MOMP is regarded as the point of no return upon which a cell is committed to undergo apoptosis ( <a href="#">Chipuk <i>et al.</i> 2006</a> ).
--------------------	--

---

<b>Western Blotting</b>	Rabbit anti BAX antibody detects a band of approximately 21 kDa in cell lysates under reducing conditions
-------------------------	---

---

<b>References</b>	<ol style="list-style-type: none"><li>Giansanti, V. &amp; Scovassi, A.I. (2008) Cell Death: A One-Way Journey to the Graveyard. <a href="#">The Open Biology Journal (1): 27-34</a>.</li><li>Chipuk, J.E. <i>et al.</i> (2006) Mitochondrial outer membrane permeabilization during apoptosis: the innocent bystander scenario. <a href="#">Cell Death Differ. 13 (8): 1396-402</a>.</li></ol>
-------------------	--

---

<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10049 available at: <a href="https://www.bio-rad-antibodies.com/SDS/AHP2716">https://www.bio-rad-antibodies.com/SDS/AHP2716</a> 10049
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)  
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)  
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)  
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M340150:190104'

**Printed on 18 Jan 2024**

---