

Datasheet: AHP2804Z

BATCH NUMBER 156478

Description:	GOAT ANTI MOUSE CXCL16:Preservative Free
Specificity:	CXCL16
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.5 ug/ml - 2.0 ug/ml
Western Blotting	▪			0.1 ug/ml - 0.2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.1 ml distilled water

Antiserum Preparation Produced from sera of goats immunized with recombinant mouse CXCL16 and purified by affinity chromatography

Buffer Solution Phosphate buffered saline

Preservative Stabilisers None present

Carrier Free Yes

Endotoxin Level < 0.1 ng/μg of protein (< 1 EU/μg)

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
Immunogen	<i>E. coli</i> derived recombinant mouse CXCL16
External Database Links	<p>UniProt: Q8BSU2 Related reagents</p> <p>Entrez Gene: 66102 Cxcl16 Related reagents</p>
Synonyms	Srpsox
Specificity	<p>Goat anti Mouse CXCL16 antibody, recognizes C-X-C motif Chemokine 16 (CXCL16). Unlike most chemokines, CXCL16 is a membrane bound ligand anchored through transmembrane domains. Although a soluble form of CXCL16 is generated through proteolytic cleavage. Transmembrane CXCL16 is mainly expressed on macrophages, dendritic cells, B cells and various tumor cells. In this form, it acts as an adhesion molecule though the binding of its receptor CXCR6 expressed on other cells including T cells. Transmembrane CXCL16 also binds gram positive and gram negative bacteria, playing a role in innate and adaptive immune responses (Shimaoka et al. 2003).</p> <p>Soluble CXCL16 has chemoattractant activity, and elicits this response through the receptor CXCR6 and plays a role in tumor progression. CXCL16 may be used as a biomarker for the evaluation of some tumor type prognoses (Richardson et al. 2015, Shibata et al. 2020, Gooden et al. 2014).</p>
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10294 available at: https://www.bio-rad-antibodies.com/SDS/AHP2804Z 10294
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To
find a

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M376518:210203'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)