

Datasheet: AHP2811Z

Description:	RABBIT ANTI MOUSE CCL11:Preservative Free
Specificity:	CCL11
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.5 ug/ml - 2.0 ug/ml
Western Blotting	▪			0.1 ug/ml - 0.2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.1 ml distilled water
Antiserum Preparation	Produced from sera of rabbits immunized with recombinant mouse CCL11 and purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 0.1 ng/μg of protein (< 1 EU/μg)
Approx. Protein	IgG concentration 1.0 mg/ml after reconstitution

Concentrations

Immunogen	<i>E. coli</i> derived recombinant mouse CCL11
-----------	--

External Database Links

UniProt:

[P48298](#)

[Related reagents](#)

Entrez Gene:

[20292](#)

Ccl11

[Related reagents](#)

Synonyms

Scya11

Specificity

Rabbit anti Mouse CCL11 antibody, recognizes C-C motif chemokine 11 (CCL11), also known as eotaxin, eosinophil chemotactic protein and small inducible cytokine A11. CCL11 is produced by a range of cell types including epithelial cells, fibroblasts, smooth muscle cells, astrocytes, chondrocytes and tissue resident macrophages. Additionally, in the central nervous system, it is expressed by choroid plexus epithelial cells, pericytes, and microglia under inflammatory conditions ([Teixeira et al. 2018](#)).

CCL11 acts as a chemoattractant for eosinophil recruitment to inflammatory sites, including during asthma attacks ([Wen & Rothenberg 2016](#)). Upon binding of CCL11 to its receptor CCR3, a series of signaling events are triggered which can result in calcium mobilization, CD11b upregulation, MAPK activation and actin polymerization. It can also result in a rapid change of shape associated with chemotaxis and granule release ([Conroy & Williams 2001](#)).

Storage

This product is shipped at ambient temperature.
Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10294 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP2811Z>
10294

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M439128:250523'

Printed on 23 May 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)