

Datasheet: AHP2813B

Description:	RABBIT ANTI SHEEP INTERLEUKIN-22:Biotin
Specificity:	IL-22
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.05 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Sheep
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antiserum to sheep interleukin-22 was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant ovine IL-22
External Database Links	UniProt:

Specificity **Rabbit anti Sheep Interleukin-22** recognizes sheep interleukin-22 (IL-22), a 21.4 kDa regulatory cytokine belonging to the IL-10 family.

Sheep IL-22 is produced in a range of cell types that include activated T cells and dendritic cells, subsets of splenic cells and Th17 cells. Ovine IL-22 modulates the inflammatory responses to bacterial pathogens, in particular in epithelial cells of the respiratory and gastrointestinal tracts.

ELISA This product may be used as a detection reagent in a sandwich ELISA together with [AHP2813](#) as the capture reagent. [PSP007](#) may be used as a standard.

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2813B>
10040

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373555:200923'

Printed on 18 Jan 2024