

Datasheet: AHP2815

BATCH NUMBER 153897

Description:	RABBIT ANTI SHEEP TNF ALPHA
Specificity:	TNF ALPHA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml - 5.0 ug/ml
Western Blotting	▪			0.1 ug/ml - 2.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Sheep

Species Cross Reactivity

Reacts with: Goat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antiserum to sheep TNF alpha was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant ovine TNF ALPHA
External Database Links	<p>UniProt: P23383 Related reagents</p> <p>Entrez Gene: 443540 TNF Related reagents</p>
Synonyms	TNFA, TNFSF2
Specificity	<p>Rabbit anti Sheep TNF alpha recognizes sheep tumor necrosis factor (TNF) alpha, a 25.5 kDa regulatory cytokine.</p> <p>TNF alpha is produced by activated macrophages, monocytes, endothelial cells, fibroblast and both T and B lymphocytes. Production of ovine TNF alpha can be induced by a variety of factors, dependant upon cell type, and includes bacterial toxins, IL-1, PDGF, IFN beta, NGF, oncostatin M and viral infections. TNF alpha functions in the regulation of immune cells, is able to induce fever, apoptotic cell death, sepsis, cachexia, inflammation and to inhibit tumorigenesis and viral replication.</p>
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with AHP2815B as the detection reagent. PSP005 may be used as a standard.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AHP2815 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
 Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
 Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M373558:200923'

Printed on 21 Mar 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)