

Datasheet: AHP2817B

BATCH NUMBER 153902

Description:	GOAT ANTI BOVINE INTERLEUKIN-8:Biotin
Specificity:	IL-8
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.01 ug/ml - 0.2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Dog

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG conjugated to Biotin - liquid

Antiserum Preparation

Antiserum to bovine interleukin-8 was raised by repeated immunization of goat with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Approx. Protein

IgG concentration 1.0 mg/ml

Concentrations

Immunogen Recombinant bovine IL-8

External Database Links

UniProt:

[P79255](#) [Related reagents](#)

Entrez Gene:

[280828](#) IL8 [Related reagents](#)

Synonyms CXCL8

Specificity **Goat anti Bovine Interleukin-8** recognizes bovine interleukin-8 (IL-8), also known as C-X-C Motif Chemokine Ligand 8 (CXCL8).

Bovine IL-8 is produced by macrophages and released from several cell types in response to an inflammatory stimulus. IL-8 is a member of the chemokine family and as chemotactic factor attracts neutrophils, basophils, and T cells but not monocytes.

ELISA This product may be used as a detection reagent in a sandwich ELISA together with [AHP2817](#) as the capture reagent. [PBP039](#) may be used as a standard.

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2817B>
10040

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373563:200923'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)