

Datasheet: AHP2823

**BATCH NUMBER 153914**

<b>Description:</b>	RABBIT ANTI DOG INTERLEUKIN-13
<b>Specificity:</b>	IL-13
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			0.1 ug/ml - 2.0 ug/ml
Electrophoresis	▪			1.0 ug/ml - 5.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Dog

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antiserum to dog interleukin-13 was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

### Approx. Protein Concentrations

IgG concentration 1.0 mg/ml

### Immunogen

Recombinant dog IL-13

**External Database  
Links**

**UniProt:**

[Q9N0W9](#)   [Related reagents](#)

**Entrez Gene:**

[442990](#) IL13   [Related reagents](#)

---

**Specificity**

**Rabbit anti Dog Interleukin-13** recognizes dog interleukin-13 (IL-13).

Dog interleukin-13 is an immunoregulatory cytokine produced by activated T lymphocytes. In dogs IL-13 is thought to have a regulatory role in allergic diseases. In in vitro experiments IL-13 was able to stimulate production of allergen-specific-IgE.

---

**ELISA**

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2823B](#) as the detection reagent. [PDP020](#) may be used as a standard.

---

**Storage**

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/AHP2823>  
10040

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)   [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...)   [RPE](#)

Sheep Anti Rabbit IgG (STAR36...)   [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

**North & South**   Tel: +1 800 265 7376

**America**   Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M373582:200923'

Printed on 12 Aug 2023