

Datasheet: AHP2823B

BATCH NUMBER 153915

Description:	RABBIT ANTI DOG INTERLEUKIN-13:Biotin
Specificity:	IL-13
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			0.05 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antiserum to dog interleukin-13 was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant dog IL-13
External Database Links	UniProt:

[Q9N0W9](#)[Related reagents](#)**Entrez Gene:**[442990](#)

IL13

[Related reagents](#)**Specificity****Rabbit anti Dog Interleukin-13** recognizes dog interleukin-13 (IL-13).

Dog interleukin-13 is an immunoregulatory cytokine produced by activated T lymphocytes. In dogs IL-13 is thought to have a regulatory role in allergic diseases. In vitro experiments IL-13 was able to stimulate production of allergen-specific-IgE.

ELISA

This product may be used as a detection reagent in a sandwich ELISA together with [AHP2823](#) as the capture reagent. [PDP020](#) may be used as a standard.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2823B>
10040

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373583:200923'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)