

Datasheet: AHP2824B

Description:	RABBIT ANTI DOG LEUKEMIA INHIBITORY FACTOR:Biotin
Specificity:	LEUKEMIA INHIBITORY FACTOR
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.01 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antiserum to dog leukemia inhibitory factor was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant dog Leukemia Inhibitory Factor
External Database Links	UniProt:

Specificity **Rabbit anti Dog Leukemia Inhibitory Factor** recognizes dog leukemia inhibitory factor (LIF).

Dog LIF is a chemokine belonging to the IL-6 cytokine family. LIF inhibits differentiation in normal stem cells, but in leukemic myeloid cells, it induces differentiation thereby preventing further growth.

ELISA This product may be used as a detection reagent in a sandwich ELISA together with [AHP2824](#) as the capture reagent. [PDP021](#) may be used as a standard.

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M373585:200923'

Printed on 29 Aug 2021