

Datasheet: AHP2829

BATCH NUMBER 153926

Description:	CHICKEN ANTI RABBIT GRO GAMMA
Specificity:	GRO GAMMA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.0 ug/ml - 5.0 ug/ml
Western Blotting	▪			0.1 ug/ml - 2.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to rabbit GRO gamma was raised by repeated immunization of chicken with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant rabbit GRO Gamma

**External Database
Links**

UniProt:

[P47854](#) [Related reagents](#)

Entrez Gene:

[100009192](#) GRO-G [Related reagents](#)

Specificity

Chicken anti Rabbit GRO Gamma (GRO Gamma) recognizes rabbit GRO gamma, also known as Chemokine (C-X-C motif) ligand 3 (CXCL3).

CXC3 belongs to the CXC chemokine family involved in host defence and immune system homeostasis. It controls migration and adhesion of monocytes to the endothelium.

ELISA

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2829B](#) as the detection reagent. [PRABP014](#) may be used as a standard.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP2829>
10040

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373595:200923'

Printed on 18 Jan 2024