

Datasheet: AHP2835B

BATCH NUMBER 153939

Description:	RABBIT ANTI CAT TNF ALPHA:Biotin
Specificity:	TNF ALPHA
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.05 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Cat

Species Cross Reactivity

Reacts with: Dog

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG conjugated to Biotin - liquid

Antiserum Preparation

Antiserum to cat TNF alpha was raised by repeated immunization of rabbit with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

Approx. Protein

IgG concentration 1.0 mg/ml

Concentrations

Immunogen Recombinant cat TNF ALPHA

External Database Links

UniProt:

[A0A0U5EMB8](#)

[Related reagents](#)

Specificity

Rabbit anti Cat TNF alpha recognizes cat tumor necrosis factor (TNF) alpha regulatory cytokine, also known as TNFSF2.

TNF alpha is produced by activated macrophages, monocytes, endothelial cells, fibroblast, and both T and B lymphocytes. TNF alpha functions in the regulation of immune cells, is able to induce fever, apoptotic cell death, sepsis, cachexia, inflammation and to inhibit tumorigenesis and viral replication.

ELISA

This product may be used as a detection reagent in a sandwich ELISA together with [AHP2835](#) as the capture reagent. [PFP005](#) may be used as a standard.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2835B>
10040

Regulatory

For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M373608:200923'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)