

Datasheet: AHP2839Z

Description:	RABBIT ANTI HUMAN INTERLEUKIN-25:Preservative Free
Specificity:	IL-25
Other names:	IL-17E
Format:	Preservative Free
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.5 ug/ml - 2.0 ug/ml
Western Blotting	▪			0.1 ug/ml - 0.2 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.1 ml distilled water
Antiserum Preparation	Produced from sera of rabbits immunized with recombinant human interleukin-17E and purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 0.1 ng/μg of protein (< 1 EU/μg)

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
Immunogen	<i>E. coli</i> derived recombinant human interleukin-25
External Database Links	<p>UniProt: Q9H293 Related reagents</p> <p>Entrez Gene: 64806 IL25 Related reagents</p>
Synonyms	IL17E
Specificity	<p>Rabbit anti Human interleukin-25 antibody, recognizes interleukin-25 (IL-25) also known as IL-17E. IL-25 is a member of the IL-17 family of cytokines and is expressed by various immune cells including mast cells, eosinophils, basophils and T helper type 2 cells. It is also expressed in non-hematopoietic cells such as epithelial and endothelial cells (Yao et al. 2016). As a pleiotropic cytokine, it stimulates stromal, innate immune and adaptive immune cells to regulate Th2 type inflammation (Monin & Gaffen 2018). IL-25 exerts its biological function through its receptor IL-17RB, although there is some evidence that it may also bind to the receptor IL-17RA (Yao et al. 2016). Through this interaction, IL-25 is able to induce expression of IL-4, IL-5, IL-13 and TSLP (McGeachy et al. 2019).</p> <p>IL-25 is upregulated in asthmatic conditions and can induce typical eosinophil inflammation and airway hyperresponsiveness (Yao et al. 2016). IL-25 levels are also associated with the severity of chronic allergic reactions (Monin & Gaffen 2018).</p>
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10294 available at: 10294: https://www.bio-rad-antibodies.com/uploads/MSDS/10294.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M376339:210127'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)