

Datasheet: AHP2979

**BATCH NUMBER 171776**

<b>Description:</b>	GOAT ANTI DOG INTERLEUKIN-1RA
<b>Specificity:</b>	IL-1RA
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5ug/ml - 5ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			1ug/ml - 10ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Dog

### Product Form

Purified IgG - liquid

### Antiserum Preparation

Antisera to canine interleukin-1RA were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

### Approx. Protein

IgG concentration 1.0 mg/ml

## Concentrations

---

**Immunogen** Recombinant canine interleukin-1RA

---

## External Database Links

### UniProt:

[Q9BEH0](#) [Related reagents](#)

### Entrez Gene:

[403660](#) IL1RN [Related reagents](#)

---

## Specificity

**Goat anti Dog interleukin-1RA antibody** recognizes canine interleukin-1 receptor antagonist, a 16.7 kDa protein belonging to the IL-1 cytokine family.

Canine IL-1RA is secreted by several cell types, including epithelial cells and adipocytes. IL-1RA functions by inhibiting the pro-inflammatory effects of both IL-1 $\alpha$  and IL-1 $\beta$ .

---

## ELISA

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2979B](#) as the detection reagent. Recombinant Dog interleukin-1RA ([PDP012](#)) may be used as a standard.

---

## Storage

This product is shipped at ambient temperature.

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

## Guarantee

12 months from date of despatch

---

## Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2979>

---

## Regulatory

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M439190:250523'

Printed on 29 Jan 2026

---