

## Datasheet: AHP2979B

<b>Description:</b>	GOAT ANTI DOG INTERLEUKIN-1RA:Biotin
<b>Specificity:</b>	IL-1RA
<b>Format:</b>	Biotin
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.01ug/ml - 0.2ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Dog
<b>Product Form</b>	Purified IgG conjugated to Biotin - liquid
<b>Antiserum Preparation</b>	Antisera to canine interleukin-1RA were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml

**Immunogen** Recombinant canine interleukin-1RA

---

**External Database Links**

**UniProt:**  
[Q9BEH0](#) [Related reagents](#)

**Entrez Gene:**  
[403660](#) IL1RN [Related reagents](#)

---

**Specificity**

**Goat anti Dog interleukin-1RA antibody** recognizes canine interleukin-1 receptor antagonist, a 16.7 kDa protein belonging to the IL-1 cytokine family.

Canine IL-1RA is secreted by several cell types, including epithelial cells and adipocytes. IL-1RA functions by inhibiting the pro-inflammatory effects of both IL-1 $\alpha$  and IL-1 $\beta$ .

---

**ELISA**

This product may be used as a detection antibody in a sandwich ELISA together with [AHP2979](#) as the capture reagent. Recombinant Dog interleukin-1RA ([PDP012](#)) may be used as a standard.

---

**Storage**

Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/AHP2979B>  
10040

---

**Regulatory**

For research purposes only

---

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334224:181127'

Printed on 25 Mar 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)