

Datasheet: AHP2990

| | |
|----------------------|-------------------------------------|
| Description: | CHICKEN ANTI RABBIT INTERLEUKIN-17A |
| Specificity: | IL-17A |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 0.1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | | | ▪ | |
| Immunohistology - Frozen | | | ▪ | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | ▪ | | | 1ug/ml - 5ug/ml |
| Immunoprecipitation | | | ▪ | |
| Western Blotting | ▪ | | | 0.1ug/ml - 2ug/ml |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Rabbit

Species Cross Reactivity

Reacts with: Bovine, Dog, Horse, Sheep, Rat, Human

Does not react with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antisera to rabbit interleukin-17A were raised by repeated immunisation of chickens with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

| | |
|---------------------------------------|---|
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) |
| Approx. Protein Concentrations | IgG concentration 1.0 mg/ml |
| Immunogen | Recombinant rabbit interleukin-17A |
| External Database Links | UniProt: G1SLF2 Related reagents |
| Specificity | Chicken anti Rabbit Interleukin-17A antibody recognizes rabbit interleukin-17A (IL-17A), a 35 kDa disulfide-linked homodimer consisting of two 153 amino acid chains, which belongs to the IL-17 family (classified as IL-17A to IL-17F), and signals through the IL-17 receptor (IL-17R/CDw217). IL-17, originally thought to have been produced primarily by activated CD4+ T cells, appears to have a range of cellular origins and induces a variety of target cells to secrete inflammatory cytokines, including IL-1β, IL-6, TNFα, IFNγ and granulocyte colony-stimulating factor, as well as synergizing with other inflammatory cytokines and agonists. |
| ELISA | This product may be used as a capture antibody in a sandwich ELISA together with AHP2990B as the detection reagent. Recombinant rabbit interleukin-17A (PRABP004) may be used as a standard. |
| Storage | Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use. |
| Guarantee | 12 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf |
| Regulatory | For research purposes only |

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M334235:181127'

Printed on 21 Mar 2022