

## Datasheet: AHP2991B

<b>Description:</b>	GOAT ANTI RABBIT INTERLEUKIN-4:Biotin
<b>Specificity:</b>	IL-4
<b>Format:</b>	Biotin
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.01ug/ml - 0.2ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Purified IgG conjugated to Biotin - liquid
<b>Antiserum Preparation</b>	Antisera to rabbit interleukin-4 were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml

**Immunogen** Recombinant rabbit interleukin-4

---

**External Database Links**

**UniProt:**  
[Q9MZR8](#) [Related reagents](#)

**Entrez Gene:**  
[100302454](#) IL4 [Related reagents](#)

---

**Specificity**

**Goat anti Rabbit interleukin-4 antibody** recognizes rabbit interleukin-4 (IL-4), a 16.6 kDa regulatory cytokine belonging to the IL-4 / IL-13 family. IL-4 is expressed by a range of cells including B cells, macrophages and monocytes. It plays an important role in the regulation of rabbit B cell responses, is involved in promoting secretion and cell surface expression of IgG1 and IgE, and with the regulation of CD23 expression on lymphocytes and monocytes.

---

**ELISA**

This product may be used as a detection antibody in a sandwich ELISA together with [AHP2991](#) as the capture reagent. Recombinant Rabbit interleukin-4 ([PRABP005](#)) may be used as a standard.

---

**Storage**

Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:  
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

---

**Regulatory**

For research purposes only

---

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M334238:181127'

Printed on 21 Sep 2021

---