

Datasheet: AHP2993

BATCH NUMBER 154784

Description:	GOAT ANTI RABBIT IP-10
Specificity:	IP-10
Other names:	CXCL10
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5ug/ml - 5ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1ug/ml - 2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to rabbit IP-10 were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein	IgG concentration 1.0 mg/ml

Concentrations

Immunogen Recombinant rabbit IP-10 (CXCL10)

External Database Links

UniProt:

[G1TBA6](#)

[Related reagents](#)

Specificity

Goat anti Rabbit IP-10 antibody recognizes rabbit (lapine) interferon gamma-inducible protein 10 (IP-10), also known as CXCL10, a 10.82 kDa protein belonging to the chemokine CXC family. IP-10 is expressed in a variety of tissues including liver, gut, spleen, kidney and thymus, in response to interferon gamma stimulation.

IP-10 acts as a chemoattractant for monocytes, NK cells and certain T-cells.

ELISA

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2993B](#) as the detection reagent. [PRABP007](#) may be used as a standard.

Further Reading

1. Müller, M. *et al.* (2010) Review: The chemokine receptor CXCR3 and its ligands CXCL9, CXCL10 and CXCL11 in neuroimmunity--a tale of conflict and conundrum. [Neuropathol Appl Neurobiol. 36: 368-87.](#)

Storage

Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2993>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M364184:200529'

Printed on 01 Mar 2024
