

## Datasheet: AHP2995B

<b>Description:</b>	CHICKEN ANTI RABBIT EXODUS-2:Biotin
<b>Specificity:</b>	EXODUS-2
<b>Other names:</b>	CCL21
<b>Format:</b>	Biotin
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.05ug/ml - 0.5ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Purified IgG conjugated to Biotin - liquid
<b>Antiserum Preparation</b>	Antisera to rabbit Exodus-2 (CCL21) were raised by repeated immunisation of chickens with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.02% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein</b>	IgG concentration 1.0 mg/ml

## Concentrations

---

**Immunogen** Recombinant rabbit CCL21 / Exodus-2

---

**External Database Links**

**UniProt:**

[B0ZYR1](#)

[Related reagents](#)

---

**Specificity**

**Chicken anti Rabbit Exodus-2 antibody** recognizes rabbit Exodus-2, also widely known as C-C motif chemokine 21 (CCL21), an 111 amino acid CC chemokine which, like CCL19 (small inducible cytokine A19), signals through the CCR7 receptor.

Exodus-2 is expressed in the HEV and T cell rich areas of lymph nodes and Peyer's patches, T cell areas of the spleen and by lymphatic endothelium of many organs; acting as a potent chemoattractant for thymocytes and T lymphocytes. Exodus-2 possesses characteristics of a homing molecule for lymphocytes to secondary lymphoid organs and has been shown to inhibit hematopoiesis.

---

**ELISA**

This product may be used as a detection antibody in a sandwich ELISA together with [AHP2995](#) as the capture reagent. Recombinant Rabbit Exodus-2 ([PRABP009](#)) may be used as a standard.

---

**Storage**

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP2995B>  
10040

---

**Regulatory**

For research purposes only

---

**North & South America** Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M420745:230706'

Printed on 12 Aug 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)