

Datasheet: AHP2996B

Description:	RABBIT ANTI BOVINE MIP-1 ALPHA:Biotin
Specificity:	MIP-1 ALPHA
Other names:	CCL3
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.05ug/ml - 0.5ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to bovine MIP-1alpha were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen Recombinant bovine MIP-1alpha

External Database Links

UniProt:

[Q8SQA6](#) [Related reagents](#)

Entrez Gene:

[282170](#) CCL3 [Related reagents](#)

Synonyms SCYA3

Specificity **Rabbit anti Bovine MIP-1 alpha antibody** recognises bovine MIP-1 alpha, otherwise known as CCL3, a 93 amino acid member of the CC chemokine family which plays a critical role in hematopoiesis and is shown to signal through the chemokine receptors CCR1, CCR4 and CCR5.

MIP-1 alpha is also involved in host defence against invading pathogens, acting as a regulator of activation and trafficking of inflammatory cells, such as macrophages, NK cells, granulocytes, dendritic cells and T-cells.

ELISA This product may be used as a capture antibody in a sandwich ELISA together with [AHP2996B](#) as the detection reagent. Recombinant bovine MIP-1 ([PBP038](#)) may be used as a standard.

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334248:181127'

Printed on 29 Aug 2021
