

Datasheet: AHP3007

Description:	RABBIT ANTI PHF6
Specificity:	PHF6
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Antiserum Preparation

Antiserum to PHF6 was raised by repeated immunization of rabbits with highly purified antigen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.02% Sodium Azide
50% Glycerol

Immunogen

Recombinant protein of human PHF6

**External Database
Links**

UniProt:

[Q8IWS0](#) [Related reagents](#)

Entrez Gene:

[84295](#) PHF6 [Related reagents](#)

Synonyms

KIAA1823

Specificity

Rabbit anti PHF6 antibody recognizes PHD finger protein 6, also known as CENP-31, BFLS, BORJ and KIAA1823.

The gene encoding PHF6 is located on the x-chromosome, and mutations cause Borjeson-Forssman-Lehmann Syndrome, a rare disease characterized by features including intellectual disability, developmental delay, and epilepsy (Mangelsdorf et al. 2009). PHF6 is a tumor suppressor, and a negative regulator of hemopoietic stem cell (HSC) self-renewal (McRae et al. 2019). Therefore, mutations that inactivate PHF6 gene function have also been implicated in hematological malignancies (Miyagi et al. 2019).

Storage

Store at -20°C
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10049 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP3007>
10049

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®680](#), [DyLight®800](#)

**North & South
America**

Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M361995:200415'

