

Datasheet: AHP300G

BATCH NUMBER 149813

Description:	RABBIT ANTI HUMAN INTERLEUKIN-15
Specificity:	IL-15
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA	▪			0.5ug/ml
Immunoprecipitation	▪			
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.15 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species

Human

Product Form

Purified IgG - lyophilised

Reconstitution

Reconstitute with 0.1 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.

Antiserum Preparation

Antisera to human IL-15 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by antigen affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Recombinant human IL-15.
External Database Links	<p>UniProt: P40933 Related reagents</p> <p>Entrez Gene: 3600 IL15 Related reagents</p>
RRID	AB_2124584
Specificity	<p>Rabbit anti Human Interleukin-15 antibody recognizes human IL-15 (Interleukin-15), a 14kDa immunomodulatory glycoprotein expressed by multiple tissues and cell types, which primarily stimulates T lymphocytes, but is also critical for the development and activation of natural killer cells (NK).</p> <p>IL-15 shares functional properties with IL-2 and binds with high affinity to the IL-15 receptor (IL-15R), consisting of the IL-15R alpha chain, and also the IL-2Rbeta and common gamma chains used by the IL-2 receptor. Affects of IL-15 have been linked with cells and tissues unrelated to the immune system, and increased expression of IL-15 has been demonstrated in Rheumatoid Arthritis (RA), ulcerative colitis (UC) Chron's disease (CD) and lymphoid malignancies, and also conditions affiliated with HIV and HTLV-1 infection.</p> <p>Rabbit anti Human Interleukin-15 antibody will neutralise IL-15 activity in bioassay systems.</p>
Instructions For Use	For functional studies do not use Azide.
References	<ol style="list-style-type: none"> 1. Kurys, G. <i>et al.</i> (2000) The long signal peptide isoform and its alternative processing direct the intracellular trafficking of interleukin-15. J. Biol. Chem. 275: 30653 - 30659. 2. Fehniger, T.A. and Caligiuri, M.A. <i>et al</i> (2001) Interleukin 15: biology and relevance to human disease. BLOOD. 97: 14-32.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.</p>

Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10162 available at:
<https://www.bio-rad-antibodies.com/SDS/AHP300G>
10162

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M372528:200706'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)