

Datasheet: AHP3013

BATCH NUMBER 161864

Description:	RABBIT ANTI SARS-CoV-2 SPIKE PROTEIN S1
Specificity:	SARS-CoV-2 SPIKE PROTEIN S1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			0.5 ug/ml
ELISA	▪			1000 ng/ml - 8 ng/ml
Western Blotting	▪			1 ug/ml
Immunocytochemistry	▪			5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Viral
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	1.0 mg/ml
Immunogen	Peptide corresponding to 16 amino acids near the amino terminus of SARS-CoV-2 (COVID-19) Spike S1 glycoprotein.

The immunogen is located within the first 50 amino acids of SARS-CoV-2 (COVID-19) Spike S1 protein.

External Database Links

UniProt:

[P0DTC2](#) [Related reagents](#)

Specificity

Rabbit anti SARS-CoV-2 spike protein S1 antibody recognizes spike glycoprotein S1, also known as S glycoprotein, peplomer protein and E2.

Human coronaviruses are pathogens that can cause severe respiratory tract infections. SARS-CoV-2 is a novel human coronavirus characterized in 2020, and causes "coronavirus disease 2019" or "COVID-19" (Xu et al. 2020). The spike glycoproteins of coronaviruses are incorporated into the viral envelope and promote entry into cells (Walls et al. 2020). The spike glycoprotein of SARS-CoV-2 is comprised of two subunits called S1 and S2. The S1 protein binds to a receptor on human cells, and the transmembrane unit S2 helps fuse the viral membrane with a cellular membrane (Hoffmann et al. 2020). Specifically, both SARS-CoV and SARS-CoV-2 spike proteins contain a receptor-binding domain (RBD) that recognize the same receptor, ACE2, on human cells (Shang et al. 2020). This antibody recognizes the spike S1 protein of SARS-CoV-2.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP3013>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M391383:211012'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)