

## Datasheet: AHP3015

**BATCH NUMBER 163660**

<b>Description:</b>	RABBIT ANTI SARS-CoV-2 SPIKE PROTEIN
<b>Specificity:</b>	SARS-CoV-2 SPIKE PROTEIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	IgG
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			0.2 ug/ml
ELISA	▪			1000 ng/ml - 0.5 ng/ml
Western Blotting	▪			1 ug/ml
Immunofluorescence	▪			1 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Viral
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.02% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	1.0 mg/ml
<b>Immunogen</b>	Peptide corresponding to 20 amino acids near the carboxy terminus of SARS-CoV-2 (COVID-19) Spike glycoprotein.

The immunogen is located within the last 50 amino acids of SARS-CoV-2 (COVID-19) Spike protein.

---

**External Database Links**

**UniProt:**

[P0DTC2](#)   [Related reagents](#)

---

**Specificity**

**Rabbit anti SARS-CoV-2 spike protein antibody** recognizes spike glycoprotein, also known as S glycoprotein, peplomer protein and E2.

Human coronaviruses are pathogens that can cause severe respiratory tract infections. SARS-CoV-2 is a novel human coronavirus characterized in 2020, and causes "coronavirus disease 2019" or "COVID-19" (Xu et al. 2020). The spike glycoproteins of coronaviruses are incorporated into the viral envelope and promote entry into cells (Walls et al. 2020). The spike glycoprotein of SARS-CoV-2 is comprised of two subunits called S1 and S2. The S1 protein binds to a receptor on human cells, and the transmembrane unit S2 helps fuse the viral membrane with a cellular membrane (Hoffmann et al. 2020). Specifically, both SARS-CoV and SARS-CoV-2 spike proteins contain a receptor-binding domain (RBD) that recognize the same receptor, ACE2, on human cells (Shang et al. 2020).

---

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/AHP3015>  
10040

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)   [FITC](#)

Goat Anti Rabbit IgG (Fc) (STAR121...)   [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (STAR35...)   [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

**North & South America**   Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**   Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**   Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M391381:211012'

**Printed on 19 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)