

Datasheet: AHP3114

BATCH NUMBER 171482

Description:	RABBIT ANTI HUMAN XRCC1
Specificity:	XRCC1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/500 - 1/1000
Immunofluorescence	▪			1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

Purified IgG - liquid

Antiserum Preparation

Antiserum to human XRCC1 was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.035% Sodium Azide (NaN₃)
30% Glycerol

Carrier Free

Yes

Immunogen

A recombinant protein corresponding to the N-terminal region of human XRCC1

External Database Links

UniProt:

Entrez Gene:

[7515](#) XRCC1 [Related reagents](#)

Specificity **Rabbit anti Human XRCC1 antibody** recognizes x-ray repair cross-complementing protein 1, also known as DNA repair protein XRCC1.

XRCC1 has been shown to play a role in various DNA repair mechanisms including, single-strand break repair, base excision repair and nucleotide excision repair ([London, 2015](#)). XRCC1 has been shown to play a role in various DNA repair mechanisms including, single-strand break repair, base excision repair and nucleotide excision repair ([Kang *et al.* 2010](#)).

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10049 available at: <https://www.bio-rad-antibodies.com/SDS/AHP3114>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M394936:220218'

Printed on 29 Jan 2026