

## Datasheet: AHP3115

**BATCH NUMBER 166733**

<b>Description:</b>	RABBIT ANTI RAT CALCITONIN RECEPTOR
<b>Specificity:</b>	CALCITONIN RECEPTOR
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			1/1,000 - 1/5,000
ELISA	▪			1/20,000
Western Blotting	▪			
Immunocytochemistry	▪			1/1,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Target Species

Rat

#### Species Cross Reactivity

Reacts with: Mouse

Reacts weakly with: Human

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - lyophilized

#### Reconstitution

Reconstitute with 100ul distilled water

#### Antiserum Preparation

Antisera to calcitonin receptor were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on Protein A.

<b>Buffer Solution</b>	10 mM Sodium Phosphate
<b>Preservative Stabilisers</b>	0.01% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml after reconstitution
<b>Immunogen</b>	Synthetic peptide derived from a C-terminal sequence of rat calcitonin receptor.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P32214</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">116506</a>    Calcrcr    <a href="#">Related reagents</a></p>
<b>Specificity</b>	<p><b>Rabbit anti rat calcitonin receptor antibody</b> recognizes rat calcitonin receptor, also known as CT-R or C1A/C1B.</p> <p>Calcitonin receptor is a seven-transmembrane domain G protein-coupled receptor that binds the protein hormone calcitonin. It is expressed in kidneys, central and peripheral nervous system, skeletal muscle stem cells and osteocytes. Investigations on calcitonin receptor's function have mainly focused on its role in regulating calcium homeostasis and osteoclast-mediated bone resorption. However, a growing body of evidence suggests that calcitonin receptor could also be involved in foetal morphogenesis and cancer development.</p> <p>Rabbit anti rat calcitonin receptor antibody binds to an intracellular epitope of the calcitonin receptor that is shared by both C1a and C1b rat calcitonin receptor isoforms. It has been shown to cross-react with mouse calcitonin receptor.</p>
<b>Histology Positive Control Tissue</b>	Kidney, cultured osteoclasts
<b>Storage</b>	<p>Prior to reconstitution store at +4°C.  After reconstitution store at -20°C.  Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #20505 available at:  <a href="https://www.bio-rad-antibodies.com/SDS/AHP3115">https://www.bio-rad-antibodies.com/SDS/AHP3115</a>  20505</p>
<b>Regulatory</b>	For research purposes only

## Related Products

## Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M415238:221220'

Printed on 19 Jan 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)