

Datasheet: AHP340G

BATCH NUMBER 157571

Description:	RABBIT ANTI HUMAN RANTES
Specificity:	RANTES
Other names:	CCL5
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA	▪			0.5 ug/ml - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1ug/ml - 0.2ug/ml
Functional Assays	▪			5.0ug/ml - 7.0ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.1 ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. NB: For functional studies do not add azide.

Antiserum Preparation Antisera to human RANTES were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by antigen affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
Immunogen	Recombinant Human RANTES
External Database Links	<p>UniProt: P13501 Related reagents</p> <p>Entrez Gene: 6352 CCL5 Related reagents</p>
Synonyms	D17S136E, SCYA5
RRID	AB_2228621
Specificity	<p>Rabbit anti Human RANTES antibody recognizes human RANTES (regulated upon activation, normal T cell expressed and secreted), otherwise known as CCL5, a 68 amino acid CC- chemokine involved in T cell activation and cell recruitment, expressed by cytotoxic (CD8+) T cells (CTLs), macrophages, keratinocytes and endothelial cells.</p> <p>RANTES acts as a chemoattractant for monocytes, CD4+/CD45RO+ memory T cells, basophils and eosinophils, thorough binding to chemokine (C-C motif) receptor 1 (CCR1), CCR3 and CCR5. Both CCR3 and CCR5 are known co-receptors for the entry of HIV viruses into target cells and the binding of RANTES to CCR3 on CTLs has been shown to enhance HIV-1-specific cytotoxicity, whilst binding to CCR5 inhibits membrane fusion and hence HIV entry.</p>
ELISA	This purified human RANTES antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with our biotinylated human RANTES antibody (AHP340B) as the detection reagent and recombinant human RANTES (PHP062) as the standard.
Western Blotting	This antibody may be used in Western Blotting under either reducing or non-reducing conditions with recombinant human RANTES (PHP062) as the positive control.
References	<ol style="list-style-type: none"> Appay, V. <i>et al.</i> (2000) RANTES activates antigen-specific cytotoxic T lymphocytes in a mitogen-like manner through cell surface aggregation. Int Immunol. 12 (8): 1173-82. Dragic, T. <i>et al.</i> (1996) HIV-1 entry into CD4+ cells is mediated by the chemokine receptor CC-CKR-5. Nature. 381 (6584): 667-73. Cheng, H. <i>et al.</i> (2016) TWEAK/Fn14 activation induces keratinocyte proliferation under psoriatic inflammation. Exp Dermatol. 25 (1): 32-7.

Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	6 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/AHP340G 10162</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	<p>Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com</p>	Worldwide	<p>Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com</p>	Europe	<p>Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com</p>
----------------------------------	--	------------------	--	---------------	--

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M372529:200706'

Printed on 28 May 2025