

Datasheet: AHP420

Description:	RABBIT ANTI LAMININ
Specificity:	LAMININ
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.25 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			1/200 - 1/1000
Immunohistology - Paraffin	▪			
ELISA	▪			>1/1000
Immunoprecipitation			▪	
Western Blotting			▪	
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Species Cross Reactivity

Reacts with: Rabbit, Chicken, Duck, Human, Dog, Pigeon

Based on sequence similarity, is expected to react with: Birds

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Whole serum - liquid

Antiserum Preparation

Antisera to laminin were raised by repeated immunisation of rabbits with highly purified antigen.

Preservative Stabilisers

0.09% Sodium Azide

Immunogen	Purified laminin from Engelbreth-Holm-Swarm (EHS) murine sarcoma.
RRID	AB_2133777
Specificity	Rabbit anti Laminin antibody recognizes both A (~400 kDa) and B (~220 kDa) chains of laminin.
Histology Positive Control Tissue	Human kidney

References

- Gordon-Weeks, P.R. *et al.* (1989) Transient expression of laminin immunoreactivity in the developing rat hippocampus. [J Neurocytol. 18 \(4\): 451-63.](#)
- Gordon-Weeks, P.R. *et al.* (1992) A study of the expression of laminin in the spinal cord of the frog during development and regeneration. [Exp Physiol. 77 \(5\): 681-92.](#)
- Abd-Elmaksoud, A. *et al.* (2009) Comparative expression of laminin and smooth muscle actin in the testis and epididymis of poultry and rabbit. [J Mol Histol. 40: 407-16.](#)
- Au, E. *et al.* (2007) SPARC from olfactory ensheathing cells stimulates Schwann cells to promote neurite outgrowth and enhances spinal cord repair. [J Neurosci. 27: 7208-21.](#)
- Mueller, S.N. *et al.* (2007) Viral targeting of fibroblastic reticular cells contributes to immunosuppression and persistence during chronic infection. [Proc Natl Acad Sci U S A. 104: 15430-5.](#)
- Falk, V. *et al.* (2002) Regulation of matrix metalloproteinases and effect of MMP-inhibition in heart transplant related reperfusion injury. [Eur J Cardiothorac Surg. 22: 53-8.](#)
- Richardson, G.D. *et al.* (2009) Dynamic expression of Syndecan-1 during hair follicle morphogenesis. [Gene Expr Patterns. 9: 454-60.](#)
- Miners, J.S. *et al.* (2008) Angiotensin-converting enzyme (ACE) levels and activity in Alzheimer's disease, and relationship of perivascular ACE-1 to cerebral amyloid angiopathy. [Neuropathol Appl Neurobiol. 34: 181-93.](#)
- Melo, F.A. *et al.* (2009) Hepatic extracellular matrix alterations in dogs naturally infected with Leishmania (Leishmania) chagasi. [Int J Exp Pathol. 90: 538-48.](#)
- Yeh, T.S. *et al.* (2014) Baculovirus-transduced, VEGF-expressing adipose-derived stem cell sheet for the treatment of myocardium infarction. [Biomaterials. 35 \(1\): 174-84.](#)
- Foditsch, E.E. *et al.* (2016) Skeletal muscle proteins: a new approach to delimitate the time since death. [Int J Legal Med. 130 \(2\): 433-40.](#)
- Kato, Y. *et al.* (2015) Targeting Antigen to Clec9A Primes Follicular Th Cell Memory Responses Capable of Robust Recall. [J Immunol. 195 \(3\): 1006-14.](#)
- Göransson, S. *et al.* (2023) Extracellular matrix stiffness-induced breast cancer cell transcriptome resembles the transition from ductal carcinoma in situ (DCIS) to invasive ductal carcinoma (IDC) [Biochem Biophys Res Comms 2 Mar. \[Epub ahead of print\]. Pre-proof.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10081 available at: https://www.bio-rad-antibodies.com/SDS/AHP420 10081
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M418578:230427'

Printed on 27 Apr 2023